



### DGP-PM Features

The DGP-PM Disposable Ground Pad features include:

- A large-area, 110 square centimeter dispersive reference electrode that is for standard pain therapy or neurosurgical RF procedures.
- Comes individually packed.
- Has a 9-foot leader wire with 4 mm shrouded safety connector.
- Requires no extra extension cable.
- Compatible with Cosman RF Generators.
- Compatible with Radionics RF Generators.

### Catalog Numbers

The DGP-PM's are sold in boxes of individually packaged DGP-PM's specified as follows:

**DGP-PM-5;** Box of 5 each DGP-PM's.

**DGP-PM-10;** Box of 10 each DGP-PM's.

**DGP-PM-25;** Box of 25 each DGP-PM's.

